



National Aeronautics and  
Space Administration




# **ARC Cell Science Validation (CS-V) Payload Overview**

**April 2017**

**POIWG #41**

**Nikita Gilkerson**

## Description of the Bioculture System



- Automated cell biology system for laboratory and International Space Station (ISS) National Laboratory research.
- Enhanced cell culture platform that provides undisturbed culture maintenance, including feedback temperature control, medical grade gas supply, perfusion nutrient delivery and removal of waste, and automated experiment manipulations.
- Programmable manipulations include: media feeds/change out, injections, fraction collections, fixation, flow rate, and temperature modification within a one-piece sterile barrier flow path.
- Cassette provides 3 levels of containment and allows Crew access to the bioculture chamber and flow path assembly for experiment initiation, refurbishment, or sample retrieval and preservation.



## CS-V Goal and Objectives

**Goal:** Validate the performance of the Bioculture System in the space flight environment on ISS.


**Objective 1: Biology.** The Bioculture System maintains viable cell cultures to the experiment requirements and preserves cells for post-flight analyses. Biology will use human Cardiomyocytes (four cassettes) and mouse MLO-Y4 osteocyte-like cells (four Cassettes) to validate culture maintenance hardware. The biology will be characterized and analyzed post-flight by examining cell morphology, quality of the cell conditioned medium, RNA quantity and quality, and gene expression.

**Objective 2: Engineering.** The Bioculture System operates as an automated incubator system for cellular and microbiological ISS research with capabilities to provide selectable environmental conditions. Engineering will have two dedicated Cassettes to conduct testing of the environmental and fluid subsystems.

**Objective 3: Operations.** The Bioculture System enables on-orbit crew operations for routine maintenance and experiment-specific tasks. Operations will evaluate crew procedures and on-orbit timelines, data handling, and ground commanding.



## CS-V Timeline



\* Best effort for no more than 15 minutes powered off during transfers  
- Asynchronous ground controls: 8 to 10 day delay for ground controls at ARC



## Operations Flight Overview

**Manifested for SpX-13 Launch (Nov 2017) and SpX-13 Return (increment 53/54)**

**Powered internal payload, operational for full mission duration**

**Payload Integration Agreement (PIA) Letter (signed) contains the following unique agreements:**

- Late load of L-24 hours or later (request for L-18 hours)
- 150 W powered on Dragon (2 locker spaces) - **Power interruption of no more than 15 minutes at any time**
- Early retrieval of continuously powered Bioculture System and cold stow samples at Recovery dock - **Live cells returning in the Bioculture System**
- Soft stowed items returned at JSC

**Pre-flight specimen and hardware processing in the KSC SSPF**

**Post-flight science recovery at ARC laboratory**

**Use of ISS facilities:**

- ExPRESS Rack
- Microgravity Science Glovebox (MSG)
- MELFI (and Cold Stowage on Dragon – GLACIER or DCB)
- Wetlab Items (gloves, wipes, absorbent pads)



## PIA (#1 through #4)

- #1: Requires early retrieval of the samples (returned in cold stowage) and the powered Bioculture System locker at Long Beach, CA for transport to the Science Lab at Ames Research Center in Moffett Field, CA
- #2: Late load of L-24 hours for powered locker. The Bioculture System must be mounted on the Dragon's hatch.
- #3: Powered locker required for: pre-launch, launch, on-orbit, return, and post-landing. Request to keep power interruptions to < 15 minutes.
- #4: On-orbit use of Microgravity Science Glovebox (MSG) for operations on cassette(s): remove samples and bioreactors and replace media bags.





## PIA (#5 through #8)

- #5: On-orbit use of Wet Lab Kit: sterile wipes, gloves, absorbent pads.
- #6: Use of 3 labs at KSC Space Station Processing Facility (SSPF): 2 for Engineering, 1 for Science.
- #7: The Bioculture System requires 150 W power for ascent and descent.
- #8: The Bioculture System requires the standard health and status data provided by SpaceX along with the following non-standard data (post-flight only): cabin temperature, cabin total pressure, cabin relative humidity, cabin ppCO<sub>2</sub> content, cabin ppO<sub>2</sub> content.




National Aeronautics and  
Space Administration



## Bioculture System

10 **Cassettes** each with  
independent environmental  
control



**Gas Supply**  
(Large and Small)

**Power and Control  
Module**

**ExPRESS Rack Locker**






## Cassette with Flow Path

Cassette




Cassette Lid



Media, Sump,  
Sample &  
Fixative Bags


Flow Path



Bioreactor



## CS-V On-Orbit Crew Ops





## Operational Support Hardware (OSH)

### Ascent



#### Soft Stow:

- Bioreactor Container
- Waste Bag
- Wipes Kit (for use in MSG)
- Tube Cutter
- O-rings (for use in MSG)
- Replacement Gas Supply Assembly
- CO2 canister

#### Cold Stow:

- Media Stowage Bags
- Sump Stowage Bags

**\*\*OpNom are still under development\*\***

### Return



#### Soft Stow:

- Waste Bag
- Wipes Kit (for use in MSG)
- Tube Cutter
- O-rings (for use in MSG)
- Original Gas Supply Assembly
- CO2 canister

#### Cold Stow: -20 °C or colder

- Bioreactor Container
- Bioreactors
- Media Stowage bags
- Sump Stowage Bags
- Sample Bags



## Cold Stowage Overview

### Launch – L-3 Turn Over

- Media Stowage Bag A
  - 4 Media Bags for MLYO4 Cells
- Media Stowage Bag B
  - 4 Media Bags for iCardiomyocyte Cells

### On-Orbit Usage

- Removed on FD 9
  - Media Stowage Bag A
- Added on FD 9
  - **Return Bag A**
    - 4 Used Sump Bags & 4 Used Media Bags
- Removed on FD 14
  - Media Stowage Bag B
- Added on FD 9
  - **Return Bag B**
    - 4 Used Sump Bags & 4 Used Media Bags
  - 1 Bioreactor Container:
    - 4 Bioreactors
    - 4 Sample Bags

### Return at Dock

- **Return Bag A**
  - 4 Used Sump Bags & 4 Used Media Bags
- **Return Bag B**
  - 4 Used Sump Bags & 4 Used Media Bags
- 1 Bioreactor Container:
  - 4 Bioreactors
  - 4 Sample Bags

\* Pending OpNom Approval





National Aeronautics and  
Space Administration



# Thank you!

**Mission Support team: Please join our splinter  
session**